

Accusorb MRI™ Product List

MWT Label	Description	Size (in)	Size (cm)	
MRI M-01	Large Square Blanket	36 x 36	91 x 91	
MRI M-02	*Arm Sleeves (pair)	n/a	n/a	
MRI M-03	Knee Straps (pair)	5 x 24+28	13 x 61+71	
MRI M-04	Poncho	n/a	91 x 91	
MRI M-05	*Small Torso / Leg Wrap	24 x 30	61 x 76	
MRI M-06	*Standard Torso Wrap	24 x 52	60 x 132	
MRI M-07	Breast Shield	14 x 34	36 x 86	
MRI M-08	Boot	n/a	n/a	
MRI M-09	Small Square Blanket	20 x 20	51 x 51	
MRI M-SE	Strap Extender Set (3)	8.25	23	
MRI M-XX	Special Order	n/a	n/a	

Notes:

- Items marked with a (*) are most versatile and recommended for first time buyers.
- Recommended starter kit is a pair of Arm Sleeves (M-02) and a Small torso/Leg wrap (M-05).
- Knee straps and Arm sleeves are sold as pairs unless otherwise requested.
- The knee straps in a pair are of different lengths - the upper one is longer for better fit.
- Special order price is for 2 dimensional product only. Any special design work will be additional.
- All sizes are approximate.